

# AmpliSens® WNV-FRT PCR kit



For Professional Use Only

## Instruction Manual

### KEY TO SYMBOLS USED

	Catalogue number		Caution
	Batch code		Contains sufficient for <n> tests
	<i>In vitro</i> diagnostic medical device		Use-by Date
	Version		Consult instructions for use
	Temperature limit		Keep away from sunlight
	Manufacturer		Negative control of amplification
	Date of manufacture		Negative control of extraction
	Authorized representative in the European Community		Positive control of Amplification
	Positive Control of Extraction		Internal control

### 1. INTENDED USE

AmpliSens® WNV-FRT PCR kit is an *in vitro* nucleic acid amplification test for qualitative detection of *West Nile virus (WNV)* RNA in the clinical material (blood plasma, serum, leukocytic fraction of blood, cerebrospinal fluid, and urine) and human autopsy material (brain, liver, spleen, and lymph node tissue), animal material (brain tissue), mosquitoes and ticks using real-time hybridization-fluorescence detection of amplified products.

**NOTE:** The results of PCR analysis are taken into account in complex diagnostics of disease.

### 2. PRINCIPLE OF PCR DETECTION

*West Nile virus* detection by the polymerase chain reaction (PCR) is based on the amplification of pathogen genome specific region using specific WNV primers. In the real-time PCR, the amplified product is detected with the use of fluorescent dyes. These dyes are linked to oligonucleotide probes, which bind specifically to the amplified product during thermocycling. The real-time monitoring of fluorescence intensities during the real-time PCR allows the detection of accumulating product without re-opening the reaction tubes after the PCR run.

AmpliSens® WNV-FRT PCR kit is a qualitative test that contains the Internal Control (Internal Control STI-87-rec (IC)). It must be used in the extraction procedure in order to control the extraction process of each individual sample and to identify possible reaction inhibition.

AmpliSens® WNV-FRT PCR kit uses "hot-start," which greatly reduces the frequency of nonspecifically primed reactions. "Hot-start" is guaranteed by separation of nucleotides and Taq-polymerase using a chemically modified polymerase (TaqF). The chemically modified polymerase (TaqF) is activated by heating at 95 °C for 15 min.

The results of amplification are registered in the following fluorescence channels:

Table 1

Channel for fluorophore	FAM	JOE
cDNA-target	Internal Control (IC) STI-87-rec cDNA	WNV cDNA
Target gene	Artificially synthesized sequence	5' non-coding target and part of the genome encoding nucleoprotein (5'-UTR – C protein)

### 3. CONTENT

AmpliSens® WNV-FRT PCR kit is produced in 1 form:

variant FRT, R-V53(RG,iQ,Mx)-CE

Variant FRT kit includes:

Reagent	Description	Volume, ml	Quantity
RT-G-mix-2	colorless clear liquid	0.015	1 tube
RT-PCR-mix-1-FRT WNV	clear liquid from colorless to light lilac colour	0.6	1 tube
RT-PCR-mix-2-FEP/FRT	colorless clear liquid	0.3	1 tube
Polymerase (TaqF)	colorless clear liquid	0.03	1 tube
TM-Revertase (MMIv)	colorless clear liquid	0.015	1 tube
Positive Control cDNA WNV / STI (C+ w/wv/sti)	colorless clear liquid	0.1	1 tube
RNA-buffer	colorless clear liquid	0.6	2 tubes
Negative Control (C-)*	colorless clear liquid	1.6	8 tubes
Positive Control WNV-rec	colorless clear liquid	0.03	5 tubes
Internal Control STI-87-rec (IC)**	colorless clear liquid	0.12	5 tubes

\* must be used in the extraction procedure as Negative Control of Extraction.

\*\* add 10 µl of Internal Control STI-87-rec (IC) during the RNA extraction procedure directly to the sample/lysis mixture (see RIBO-zol-C K2-13-50-CE, RIBO-prep K2-9-Et-50-CE, MAGNO-sorb K2-16-1000-CE protocols).

Variant FRT is intended for 60 reactions (including controls).

### 4. ADDITIONAL REQUIREMENTS

- 0.15 M NaCl or PBS buffer solution (137 mM sodium chloride; 2.7 mM potassium chloride; 10 mM sodium monophosphate; 2 mM potassium diphosphate; pH=7,5±0,2).
- Homogenizer for pretreatment of autopsy material, mosquitos and ticks.
- Stainless steel beads (5 and 7 mm diameter).
- RNA extraction kit.
- Disposable powder-free gloves and a laboratory coat.
- Pipettes (adjustable).
- Sterile pipette tips with aerosol filters (up to 200 µl).
- Tube racks.
- Vortex mixer.
- Desktop centrifuge with a rotor for 2-ml reaction tubes.
- PCR box.
- Real-time instruments (for example, Rotor-Gene 3000/6000 (Corbett Research, Australia); Rotor-Gene Q (QIAGEN, Germany); iCycler iQ5 (Bio-Rad, USA), Mx3000P (Stratagene, USA)).
- Disposable polypropylene tubes:
  - a) screwed or tightly closed 1.5-ml tubes
  - b) thin-walled 0.2-ml PCR tubes with domed caps if a plate-type instrument is used;
  - c) thin-walled 0.2-ml PCR tubes with flat caps if a rotor-type instrument is used.
- Refrigerator for 2–8 °C.
- Deep-freezer at the temperature from minus 24 to minus 16 °C.
- Reservoir for used tips.

### 5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol filters and use a new tip for every procedure.
- Store all extracted positive material (specimens, controls and amplicons) away from all other reagents and add it to the reaction mix in a distantly separated facility.
- Thaw all components thoroughly at room temperature before starting an assay.
- When thawed, mix the components and centrifuge briefly.
- Use disposable protective gloves and laboratory cloths, and protect eyes while samples and reagents handling. Thoroughly wash hands afterwards.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all specimens and unused reagents in accordance with local regulations.
- Samples should be considered potentially infectious and handled in biological cabinet in compliance with appropriate biosafety practices.
- Clean and disinfect all samples or reagents spills using a disinfectant, such as 0.5 % sodium hypochlorite or another suitable disinfectant.
- Avoid inhalation of vapors, samples and reagents contact with the skin, eyes, and mucous membranes. Harmful if swallowed. If these solutions come into contact, rinse the injured area immediately with water and seek medical advice if necessary.
- Safety Data Sheets (SDS) are available on request.
- Use of this product should be limited to personnel trained in DNA amplification techniques.
- Workflow in the laboratory must be one-directional, beginning in the Extraction Area and moving to the Amplification and Detection Area. Do not return samples, equipment and reagents in the area where the previous step was performed.

Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

## 6. SAMPLING AND HANDLING

**NOTE:** Obtaining samples of biological materials for PCR-analysis, transportation, and storage are described in the manufacturer's handbook [1]. It is recommended that this handbook is read before starting work.

**AmpliSens® WNV-FRT** PCR kit is intended for analysis of the RNA extracted with RNA extraction kits from the clinical material:

– **blood plasma, blood serum, and CSF**

Take a whole blood specimen in the morning after overnight fasting in a tube with 6 % EDTA solution in the ratio 1:20. Invert the closed tube several times. To collect plasma, centrifuge the tube at 1600 g for 20 min. Obtain blood serum by standard methods. Take 200 µl of clinical material for the test (if RIBO-prep kit is used for extraction) or 1 ml (if MAGNO-sorb is used for extraction).

– **leukocytic fraction of blood**

To obtain leukocytic fraction of blood (the clinical material is recommended to be used on the second week of the disease), transfer 1.5 ml of the blood with EDTA to an Eppendorf tube and centrifuge at 400 g for 10 min. Then collect 500–600 µl of plasma and centrifuge at 7,000 g for 10 min. Use the cell pellet and 200 µl of the supernatant (if extracting with RIBO-prep or combined extraction methods (RIBO-zol-C and RIBO-sorb kits or RIBO-zol-C and RIBO-prep kits) for subsequent RNA extraction.

– **internal organs of animals and autopsy material**

Homogenize internal organs of animals and autopsy material in a porcelain mortar with a pestle and prepare 10 % suspension with sterile saline solution or phosphate buffer.

If a TissueLyser LT automated homogenizer is used, set the following parameters:

- the volume of PBS buffer or 0.15 M NaCl solution for homogenizing is defined by the volume of homogenized tissue. The tissue–buffer ratio is 1 : 9; therefore, a 10 % suspension is prepared;
- the total sample volume should not exceed 1 ml for 1.5-ml tubes;
- homogenization parameters
  - for brain tissues: bead size, 5 mm; frequency, 50 Hz/s; homogenization time, 2–3 min;
  - for liver, spleen, and lymph node tissues: bead size, 7 mm; frequency, 50 Hz/s; homogenization time, 10 min;
  - for lymph nodes: bead size, 5 mm; frequency, 50 Hz/s; homogenization time, 5 min.

Take 30 µl of the suspension for RNA extraction if extracting with RIBO-prep or by a combined extraction method (RIBO-zol-C and RIBO-sorb kits) and 100 µl of the suspension if extracting by a combined extraction method (RIBO-zol-C and RIBO-prep kits).

– **urine**

Collect a urine sample to a clean vessel. If the sample will not be analyzed within 1 day after taking, transfer the urine to a 30-ml centrifuge tube or an Eppendorf tube, add glycerol (10 % of the sample volume), stir, and freeze at the temperature minus 20 °C for storage up to 1 week or at the temperature minus 70 °C for storage for a long time.

If a cooling centrifuge (temperature, 4 °C; tube volume, 30 ml; acceleration, 8,000 g) is used, follow the pretreatment algorithm described below:

Centrifuge the sample at 8,000–9,000 g for 10 min and discard the supernatant. Transfer the pellet and 1 ml of the supernatant to an Eppendorf tube. Centrifuge the sample at 8,000 g for 10 min. Discard 900 µl of the supernatant. Use the pellet with the remaining 100 µl of the supernatant for RNA extraction. In the case of a high content of salts in the urine sample, transfer 100 µl of the supernatant to an Eppendorf tube and use it for RNA extraction.

If a cooling centrifuge for 30-ml tubes and 8,000 g is absent, centrifuge 1-ml sample as described below. Extract RNA from the pellet with 100 µl of the supernatant.

– **mosquito suspension**

To prepare mosquitoes suspension, use a porcelain mortar and pestle. If a TissueLyser LT automated homogenizer is used, set the following homogenizing parameters: bead size, 5 mm; frequency, 50 Hz/s; homogenization time, 5 min; buffer volume, 700 µl (a pooled sample of 25 mosquitoes) or 1,000–1500 µl (a pooled sample of 50 mosquitoes).

Prepare mosquito pools (up to 50 mosquitoes). Homogenize gnats in sterile saline solution or phosphate buffer calculating 30 µl of the solution per 1 mosquito. Centrifuge the samples at 10,000 g for 1 minute. Collect 100 µl or the supernatant for RNA extraction.

– **tick suspension**

Form tick pools: pool of hungry ticks consists of 5–7 specimens; pool of half-full ticks consists of 2–3 specimens; pool of blood-filled ticks consists of 1 tick.

If a TissueLyser LT automated homogenizer is used, set the following homogenizing parameters for the *Hyalomma* genus ticks: bead size, 7 mm; frequency, 50 Hz/s; homogenization time, 10–12 min; buffer volume, 700 µl (hungry tick) or 1,000–1500 µl (blood-filled tick and tick pools).

To prepare tick suspension a sterile porcelain mortar and pestle can be used as well. Prior to homogenization, blood-filled ticks should be pierced to let blood out. Homogenize ticks in 700 µl (if the sample consists of 1 hungry tick) or 1–1.5 ml (if a blood-filled tick or tick pools are homogenized) of 0.15 M sodium chloride solution. Add the solution in small portions while homogenizing. Centrifuge the prepared suspension at 10,000 rpm for 1 min and use 100 µl of the supernatant for RNA extraction.

The samples are to be stored at 2–8 °C for 1 day or at the temperature not more than minus 16 °C for 1 week.

Internal organ tissues and mosquitoes are to be stored at the temperature not more than minus 16 °C for 1 week or at the temperature minus 70 °C for a long time.

Ticks are to be stored alive at most 1 month or at the temperature not more than minus 16 °C for 1 week, or at the temperature minus 70 °C for a long time.

## 7. WORKING CONDITIONS

**AmpliSens® WNV-FRT** PCR kit should be used at 18–25 °C.

## 8. PROTOCOL

### 8.1. RNA extraction

It's recommended that the following nucleic acid extraction kits are used:

Material	Extraction kit
<ul style="list-style-type: none"> <li>leukocytic fraction of blood</li> <li>suspension of internal organs</li> <li>suspension of mosquitoes</li> <li>suspension of ticks</li> <li>urine sediment (including salts)</li> </ul>	<b>RIBO-zol-C</b> <a href="#">REF</a> K2-13-50-CE (for the stage I of RNA extraction) and <b>RIBO-prep</b> <a href="#">REF</a> K2-9-Et-50-CE or <b>RIBO-sorb</b> <a href="#">REF</a> K2-1-Et-50-CE (for the stage II of RNA extraction)
<ul style="list-style-type: none"> <li>blood plasma</li> <li>blood serum</li> <li>leukocytic fraction of blood</li> <li>cerebrospinal fluid (CSF)</li> <li>internal organ homogenates</li> <li>mosquitoes homogenates</li> <li>urine sediment (without salts)</li> </ul>	<b>RIBO-prep</b> <a href="#">REF</a> K2-9-Et-50-CE
<ul style="list-style-type: none"> <li>blood plasma</li> <li>blood serum</li> <li>CSF</li> </ul>	<b>MAGNO-sorb</b> <a href="#">REF</a> K2-16-1000-CE

### 8.1.1 RIBO-prep

**RNA is extracted from blood plasma, blood serum, leukocytic fraction of blood, CSF, urine (without salt sediments), homogenates of internal organ tissues, and mosquito homogenates.**

Extract RNA according to the manufacturer's protocol taking into account next additions and improvements:

The volume of the tissue homogenate sample is 30 µl. The volume of the suspension of mosquitoes is 100 µl. The volume of the blood plasma, blood serum, CSF is 100 µl.

If extracting from leukocytic fraction of blood and urine sediment: add 300 µl of Solution for Lysis and 10 µl of Internal Control STI-87-rec (IC) directly to the tubes with the samples.

**NOTE**

Add 300 µl of Solution for Lysis and 10 µl of Internal Control STI-87-rec (IC) to the tube labeled C– (Negative Control of Extraction); Add 10 µl of Positive Control *WNV*-rec, 10 µl of Internal Control STI-87-rec (IC), and 300 µl of Solution for Lysis to the tube labeled PCE (Positive Control of Extraction).

Centrifuge the tubes under the following conditions:

- after Solution for Precipitation addition, at 10,000 g for 5 min,
- after washing with Washing Solution 3, at 10,000 g for 2 min,
- after washing with Washing Solution 4, at 10,000 g for 2 min,
- after RNA-buffer addition, at 10,000 g for 1 min.

### 8.1.2 RIBO-zol-C and RIBO-prep

**RNA is extracted from the leukocytic fraction of blood, urine sediments, suspensions of internal organ tissues, mosquitoes, and ticks.**

I stage. RIBO-zol-C

Extract RNA according to the manufacturer's protocol taking into account next additions and improvements:

The volume of the suspensions of internal organ tissues, mosquitoes, and ticks is 100 µl.

If extracting from leukocytic fraction of blood and urine sediment: add 300 µl of Solution D and 10 µl of Internal Control STI-87-rec (IC) directly to the tubes with the samples.

Add 100 µl of Negative Control (C–) and 10 µl of Internal Control STI-87-rec (IC) to the tube labeled C– (Negative Control of Extraction);

**NOTE**

Add 10 µl of Positive Control *WNV*-rec, 10 µl of Internal Control STI-87-rec (IC), and 90 µl of Negative Control (C–) to the tube labeled PCE (Positive Control of Extraction).

Centrifuge the tubes under the following conditions:

- do not centrifuge the tubes after mixing the Solution D and Internal Control STI-87-rec (IC) with the test or control samples and incubation,
- after Solution E and solution A addition centrifuge at 1,500 g for 5 s,
- after solution B addition and placing the tubes on ice (0–4 °C) for 5 min, at 10,000 g for 10 min.

The volume of the top phase collected at the final step is 200 µl.

II stage. RIBO-prep

Extract RNA according to the manufacturer's protocol taking into account next additions and improvements:

Use 200 µl of top phase obtained at I stage of RNA extraction as a test sample. Do not add Internal Control STI-87-rec (IC).

**NOTE**

Centrifuge the tubes under the following conditions:

- after Solution for Precipitation addition, at 10,000 g for 5 min,
- after washing with Washing Solution 3 and Washing Solution 4, at 10,000 g for 2 min,
- after RNA-buffer addition and incubation, at 10,000 g for 1 min.

### 8.1.3 RIBO-zol-C and RIBO-sorb

**RNA is extracted from the leukocytic fraction of blood, urine sediments, suspensions of internal organ tissues, mosquitoes, and ticks.**

I stage. RIBO-zol-C

Extract RNA according to the manufacturer's protocol taking into account next additions and improvements:

The volume of the tissue suspension sample is 30 µl. The volume of the suspension of mosquitoes and ticks is 100 µl.

If extracting from leukocytic fraction of blood and urine sediment: add 300 µl of Solution D and 10 µl of Internal Control STI-87-rec (IC) directly to the tubes with the samples.

Add 100 µl of Negative Control (C–) and 10 µl of Internal Control STI-87-rec (IC) to the tube labeled C– (Negative Control of Extraction);

Add 10 µl of Positive Control *WNV*-rec, 10 µl of Internal Control STI-87-rec (IC), and 90 µl of Negative Control (C–) to the tube labeled PCE (Positive Control of Extraction).

**NOTE**

Centrifuge the tubes under the following conditions:

- after mixing the Solution D and Internal Control STI-87-rec (IC) with the test or control samples and incubation, at 1,500 g for 5 s,
- after Solution E and solution A addition centrifuge at 1,500 g for 5 s,
- after solution B addition and placing the tubes on ice (0–4 °C) for 5 min, at 10,000 g for 10 min.

The volume of the top phase collected at the final step should be:

- 400 µl if extracting from urine sediments, homogenates of internal organs, mosquitoes, or ticks,
- 300 µl if extracting from blood sediments,
- 450 µl if extracting from Negative or Positive Controls.

II stage. RIBO-sorb

Extract RNA according to the manufacturer's protocol taking into account next additions and improvements:

Use 300, 400 or 450 µl of top phase obtained at I stage of RNA extraction as a test sample. Do not add Internal Control STI-87-rec (IC).

Use 400 µl of Lysis Solution instead of 450 µl

Centrifuge the tubes under the following conditions:

- after sorbent addition, at 1,500 g for 30 s,
- after washing with Washing Solution 1, at 1,500 g for 30 s,
- after washing with Washing Solution 3, at 5,000 g for 45 s,
- after washing with Washing Solution 4, at 6,000 g for 1 min.

After washing with Washing Solution 4 and centrifugation, incubate all tubes with open caps for 10 min at 56 °C.

After adding RNA-buffer, incubate the samples at 56 °C for 5 min and vortex every other minute.

Store an unsealed tube with RNA-buffer at the temperature not more than minus 16 °C.

### 8.1.4 MAGNO-sorb

**NOTE** Extract the RNA according to the manufacturer's protocol.

## 8.2. Preparing reverse transcription and PCR

The total reaction volume is 25 µl, the volume of RNA sample is 10 µl.

### 8.2.1 Preparing tubes for PCR

The type of tubes depends on the PCR instrument used for analysis. Use disposable filter tips for adding reagents, RNA and control samples into tubes.

1. Prepare the required number of tubes for PCR.
2. Prepare the reaction mixture for the required number of reactions. To do this, mix in a new tube RT-PCR-mix-1-FRT *WNV*, RT-PCR-mix-2-FEP/FRT, polymerase (TaqF), RT-G-mix-2, and TM-Revertase (MMV) calculating per each reaction:
  - 10 µl of RT-PCR-mix-1-FRT *WNV*;
  - 5 µl of RT-PCR-mix-2-FEP/FRT;
  - 0.5 µl of polymerase (TaqF);
  - 0.25 µl of TM-Revertase (MMV);
  - 0.25 µl of RT-G-mix-2;

Take into account that each run includes at least four control points: Positive and Negative Controls of extraction (PCE, C–) as well as Positive and Negative Controls of RT-PCR (C+*WNV/STI*, NCA).

3. Transfer 15 µl of the prepared reaction mixture to each tube.

**NOTE** Do not store the prepared mixture.

4. Using filter tips add 10 µl of RNA samples obtained at the RNA extraction stage into prepared tubes. Carefully mix by pipette.

5. Carry out the control amplification reactions:

**NCA** – Add 10 µl of RNA-buffer to the tube labeled NCA (Negative Control of Amplification).

**C+*WNV/STI*** – Add 10 µl of Positive Control cDNA *WNV/STI* to the tube labeled C+*WNV/STI* (Positive Control of Amplification).

**C–** – Add 10 µl of the sample extracted from the Negative Control (C–) reagent to the tube labeled C– (Negative control of Extraction).

**PCE** – Add 10 µl of the sample extracted from the Positive Control *WNV-rec* reagent to the tube labeled PCE (Positive control of Extraction)

**NOTE** Amplification should immediately follow after compounding of the reaction mix with RNA-samples and controls.

### 8.2.2 Reverse transcription and amplification

1. Create a temperature profile on your instrument as follows:

Table 2

Step	Rotor-type instruments <sup>1</sup>			Plate-type instruments <sup>2</sup>		
	Temperature, °C	Time	Cycles	Temperature, °C	Time	Cycles
1	50	30 min	1	50	30 min	1
2	95	15 min	1	95	15 min	1
3	95	5 s	5	95	5 s	5
	56	25 s		56	30 s	
	72	15 s		72	15 s	
4	95	5 s	40	95	5 s	40
	56	25 s		56	30 s	
		fluorescent signal detection			fluorescent signal detection	
		72			15 s	

Fluorescent signal is detected in the channels for the FAM and JOE fluorophores.

2. Adjust the fluorescence channel sensitivity according to the *Important Product Information Bulletin* and Guidelines [2].

3. Insert tubes into the reaction module of the device.

4. Run the amplification program with fluorescence detection.

5. Analyze results after the amplification program is completed.

## 9. DATA ANALYSIS

Analysis of results is performed by the software of the real-time PCR instrument used by measuring fluorescence signal accumulation in two channels:

- The signal of the IC cDNA amplification product is detected in the channel for the FAM fluorophore.
- The signal of the *WNV* cDNA amplification product is detected in the channel for the JOE fluorophore.

Results are interpreted by the crossing (or not-crossing) the fluorescence curve with the threshold line set at the specific level that corresponds to the presence (or absence) of a Ct value of the cDNA sample in the corresponding column of the results grid.

The results are interpreted according to the Table 3.

Table 3

Correspondence of targets and channels		
PCR-mix-1	Detection in the channel for the fluorophore	
	FAM	JOE
RT-PCR-mix-1-FRT <i>WNV</i>	IC	<i>WNV</i>

Principle of interpretation is the following:

- *WNV* cDNA is **detected** in a sample if the Ct value determined in the results grid in the channel for the JOE fluorophore is less than the boundary Ct value specified in the *Important Product Information Bulletin*. Moreover, the fluorescence curve of the sample should cross the threshold line in the area of typical exponential growth of fluorescence.
- *WNV* cDNA is **not detected** in a sample if the Ct value determined in the result grid in the channel for the FAM fluorophore is less than the specified boundary value, whereas the Ct value in the channel for the JOE fluorophore is not determined or greater than the specified value in the *Important Product Information Bulletin*.
- The result is **invalid** if the Ct value is not determined (absent) in the channel for JOE fluorophore, whereas the Ct value in the channel for the FAM fluorophore is not determined (absent) or greater than the specified boundary Ct value. In such cases, the PCR analysis should be repeated starting from the RNA extraction stage.

**NOTE** Boundary Ct values are specified in the *Important Product Information Bulletin* enclosed to the PCR kit. See also Guidelines [2].

The result of the analysis is considered reliable only if the results obtained for the Positive and Negative Controls of amplification as well as for the Positive and Negative Control of extraction are correct (see Table 4).

Table 4

Results for controls			
Control	Stage for control	Ct value in the channel for the fluorophore	
		JOE	FAM
C–	RNA extraction	Absent	<boundary value
PCE	RNA extraction	<boundary value	<boundary value
NCA	PCR	Absent	Absent
C+ <i>WNV/STI</i>	PCR	<boundary value	<boundary value

## 10. TROUBLESHOOTING

Results of the analysis are not taken into account in the following cases:

1. If the Ct value determined for the Positive Control of Amplification (C+) in the channel for the JOE fluorophore is greater than the boundary Ct value or absent, the amplification should be repeated for all samples in which the specific cDNA was not detected.
2. If the Ct value determined for the Positive Control of Extraction (PCE) in the channel for the JOE fluorophore is greater than the boundary Ct value or absent, the extraction should be repeated for all samples in which the specific cDNA was not detected.
3. If the Ct value is determined for the Negative Control of Extraction (C–) in the channel for the JOE fluorophore, the PCR-analysis should be repeated for all samples in which cDNA was detected in the channel for the JOE fluorophore.
4. If the Ct value is determined for the Negative Control of Amplification (NCA) in the channels for the FAM and JOE fluorophores, the amplification should be repeated for all samples in which cDNA was detected in the channel for the JOE fluorophore, with simultaneously running NCA in triplicate.

If you have any further questions or if encounter problems, please contact our Authorized representative in the European Community.

## 11. TRANSPORTATION

AmpliSens® *WNV-FRT* PCR kit should be transported at 2–8 °C for no longer than 5 days.

## 12. STABILITY AND STORAGE

All components of the AmpliSens® *WNV-FRT* PCR kit (except for RT-G-mix-2, RT-PCR-mix-1-FRT *WNV*, RT-PCR-mix-2-FEP/FRT, polymerase (TaqF), and TM-Revertase (MMV)) are to be stored at 2–8 °C when not in use. All components of the AmpliSens® *WNV-FRT* PCR kit are stable until the expiration date on the label. The shelf life of reagents before and after the first use is the same, unless otherwise stated.

**NOTE:** RT-G-mix-2, RT-PCR-mix-1-FRT *WNV*, RT-PCR-mix-2-FEP/FRT, polymerase (TaqF), and TM-Revertase (MMV) are to be stored at temperature from minus 24 to minus 16 °C when not in use.

**NOTE:** RT-PCR-mix-1-FRT *WNV* is to be kept away from light.

## 13. SPECIFICATIONS

### 13.1. Analytical sensitivity

The analytical sensitivity of AmpliSens® *WNV-FRT* PCR kit is specified in the table below.

Type of biological material (test sample volume)	RNA extraction kit	PCR kit	Analytical sensitivity, copies/ml	Pretreatment of biological material
blood serum (200 µl), CSF (200 µl), leukocytic fraction of blood (200 µl), 10 % brain tissue homogenate (30 µl), mosquitoes (100 µl)	RIBO-prep	variant FRT	5 x 10 <sup>3</sup>	Indicated sensitivity can be reached only if the specified pretreatment instructions are followed and the specified specimen volume is used
leukocytic fraction of blood (200 µl), 10 % brain tissue homogenate (30 µl), mosquitoes (100 µl)	RIBO-zol-C RIBO-prep	variant FRT	5 x 10 <sup>3</sup>	
leukocytic fraction of blood (200 µl), 10 % brain tissue homogenate (30 µl), mosquitoes (100 µl)	RIBO-zol-C RIBO-sorb	variant FRT	5 x 10 <sup>3</sup>	
blood plasma, blood serum, CSF (1 ml for all)	MAGNO-sorb	variant FRT	5 x 10 <sup>2</sup>	

### 13.2. Analytical specificity

The analytical specificity of AmpliSens® *WNV-FRT* PCR kit is ensured by selection of specific primers and probes as well as stringent reaction conditions. The primers and probes have been checked for possible homologues to all sequences deposited in gene banks by sequence comparison analysis. The analytical specificity was assayed using the following microorganisms:

- flaviviruses (*tick-borne encephalitis virus*, *Langat*, *Powassan*, *Japanese encephalitis*, and *Omsk hemorrhagic fever viruses*);
- herpes viruses (I and II types, *CMV*, *EBV*, *VZV*, IV type), *enteroviruses* (*ECHO*, *Coxsackie*);
- rickettsiae of spotted fever group (*Rickettsia conorii* spp. *caspia* and *R.heilongjiangensis*; *Coxiella burnetii*; and *Bartonella henselae* and *B. quintana*);
- spirochaetes (*Borrelia miyamotoi*; *Treponema pallidum*; *Leptospira interrogans*, *L.kirschneri*; and *L.borgpetersenii*);

No false-positive results were observed during examination of RNA/DNA of the above-mentioned organisms, human DNA, DNA of birds, DNA of ticks and mosquitoes, DNA of rodents.

The clinical specificity of AmpliSens® *WNV-FRT* PCR kit was confirmed in laboratory clinical trials.

<sup>1</sup> For example, Rotor-Gene 3000/6000 (Corbett Research, Australia), Rotor-Gene Q (Qiagen, Germany), or equivalent.

<sup>2</sup> For example, iQ5 (Bio-Rad, USA), Mx3000P (Stratagene, USA), or equivalent.

## 14. REFERENCES

1. Handbook "Sampling, Transportation, and Storage of Clinical Material for PCR diagnostics", developed by Federal Budget Institute of Science "Central Research Institute for Epidemiology" of Federal Service for Surveillance on Consumers' Rights Protection and Human Well-Being.
2. Guidelines to the **AmpliSens® WNV-FRT** PCR kit for qualitative detection of *West Nile virus* RNA in the biological material by the polymerase chain reaction (PCR) with real-time hybridization-fluorescence detection developed by Federal Budget Institute of Science "Central Research Institute for Epidemiology".

## 15. QUALITY CONTROL

In accordance with Federal Budget Institute of Science "Central Research Institute for Epidemiology" ISO 13485-Certified Quality Management System, each lot of **AmpliSens® WNV-FRT** PCR kit has been tested against predetermined specifications to ensure consistent product quality.

List of Changes Made in the Instruction Manual

VER	Location of changes	Essence of changes
13.02.11	Cover page	The phrase "For Professional Use Only" was added
	Content	New sections "Working Conditions" and "Transportation" were added
		The "Explanation of Symbols" section was renamed to "Key to Symbols Used"
	Stability and Storage	The information about the shelf life of open reagents was added
	Key to Symbols Used	The explanation of symbols was corrected
	Through the text	The reagent «Positive Control cDNA WNV» was changed into «Positive Control cDNA WNV / STI», the abbreviation «C+ <sub>WNV</sub> » was changed into «C+ <sub>WNV / STI</sub> »
	Content	The volume of Negative Control (C-) was changed from 1.2 into 1.6 ml
	Protocol, item 8.1.3 RIBO-prep	The material «brains tissue homogenates» was changed into «homogenates of internal organ tissues »
	Protocol, item 8.2.2	In the table the numeration of steps was corrected
	Protocol, item 8.1	RIBO-sorb reagent kit was deleted, material was corrected
	Sampling and handling	Single items for blood serum and CSF were deleted
		Item for blood plasma, blood serum, and CSF was corrected
		At the leukocytic fraction of blood, internal organs of animals and autopsy material items the used extraction kits were corrected
Data analysis, result interpretation	The storage conditions of internal organ tissues were added	
Specification, sensitivity	Table 2 was added	
	In the table «clinical material» was changed into «biological material»	
07.12.11 VV	Data analysis, Interpretation of results	Analytical sensitivity for RIBO-prep and RIBO-zol-C kits was added
		For RIBO-sorb and RIBO-zol-C the biological material was changed. Blood serum and CSF were deleted, the volume of leukocytic fraction of blood was changed from 100 to 200 µl.
09.06.16 ME	Text	The principle of interpretation of results was edited
	4. Additional requirements	Corrections according to the template
	8.1. RNA extraction	Additional requirements for pretreatment was added
	8.2. Preparing reverse transcription and PCR	MAGNO-sorb nucleic acid extraction kit is recommended. The section is rewritten
18.03.19 EM	3. Content	Information about controls of extraction was added
		The reference to Guidelines was added
19.05.20 MM	Through the text	The colour of the reagent was specified
	2. Principle of PCR detection	The text formatting was changed
12.03.21 MM	Footer	The table with targets was added.
		The phrase "Not for use in the Russian Federation" was added
—	—	The name, address and contact information for Authorized representative in the European Community was changed

**AmpliSens®**



Ecoli Dx, s.r.o., Purkyňova 74/2  
110 00 Praha 1, Czech Republic  
Tel.: +420 325 209 912  
Cell: +420 739 802 523



Federal Budget Institute of  
Science "Central Research  
Institute for Epidemiology"  
3A Novogireevskaya Street  
Moscow 111123 Russia