



## FavorPrep™ 96-well Viral DNA/ RNA Extraction Kit

(For Research Use Only)

- For 96-well high-throughput extraction of viral nucleic acid from cell free samples such as serum, plasma, body fluids and the supernatants of cell cultures

### Kit Contents:

Cat. No.: (Q'ty)	FAVRE 96001 (1 plate)	FAVRE 96002 (2 plates)	FAVRE 96004 (4 plates)
VNE Buffer	60 ml	120 ml	120 ml x 2
AD Buffer <sup>a</sup> (concentrate)	5 ml	10 ml	10 ml x 2
Wash Buffer 1 <sup>*</sup> (concentrate)	55 ml	110 ml	110 ml x 2
Wash Buffer 2 <sup>■</sup> (concentrate)	25 ml	50 ml	50 ml x 2
RNase-Free Water	15 ml	30 ml	30 ml x 2
Filter Plate (96-Well nucleic acid binding plate)	1 plate	2 plates	4 plates
Collection Plate (96-Well 2 ml Plate)	3 plates	6 plates	12 plates
Elution Plate (96-Well PCR plate)	1 plate	2 plates	4 plates
Adhesive Film	2 pcs	4 pcs	8 plates

### Preparation of working buffers

Add RNase-free ethanol (96~100%) to AD Buffer, Wash Buffer 1 and Wash Buffer 2 when first use.

	FAVRE 96001	FAVRE 96002	FAVRK 96004
<sup>a</sup> Ethanol volume for AD Buffer	40 ml	80 ml	
<sup>*</sup> Ethanol volume for Wash Buffer 1	10 ml	20 ml	
<sup>■</sup> Ethanol volume for Wash Buffer 2	100 ml	200 ml	

### Related products can be ordered from Favorgen

	Cat. No:	Description:
Vacuum manifold	Wel-Vac 200	Size: 23.2x12.4x10.2 cm; material: anodized aluminum Unique designed column adaptor board and luer connector make Wel-Vac 200 highly compatible with various kinds of centrifuge tube and 96-well plate.
Oil-less vacuum pump	FAPMP 110/220	FAPMP 110: 110V, 60Hz, Max. vacuum -26.8 inches Hg (-680 mm Hg) FAPMP 220: 220V, 50Hz, Max. vacuum -26.8 inches Hg (-680 mm Hg)

### Quality control

The quality of 96-Well Viral DNA/ RNA Extraction Kit is tested on a lot-to-lot basis. The purified nucleic acid is checked by real-time PCR and capillary electrophoresis,

### Specification:

Principle: Filter Plate (silica membrane)

Sample size: up to 200 µl of serum, plasma, body fluids and the supernatant of cell cultures

Processing: centrifugation protocol or vacuum & centrifugation protocol

Operation time: < within 1 hr/ 96 preparations

RNA Binding capacity: up to 75 µg/ well

Elution volume: 50 ~ 75 µl

### Reagent to be provided by user

96 ~100 % RNase - free ethanol

### Important notes:

1. Make sure everything is RNase-free when handling this kit.
2. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
3. Add RNase-free ethanol (96~100%) to AD Buffer, Wash Buffer 1 and Wash Buffer 2 when first use.
4. Equipments required:
  - For centrifugation protocol: A centrifuge is required, capable of 5,600 ~ 6,000 X g, with a swing -bucket rotor and the adaptor for 96-well plates.
  - For vacuum protocol: A vacuum manifold for 96-well plate and a vacuum source reached to 15 inches Hg are required.

**(Alternative):** If using centrifugation for Elution Step (STEP 8), a centrifuge equipment is required, capable of 5,600 ~ 6,000 X g, with a swing -bucket rotor and the adaptor for 96-well plate.

### Brief procedure:

**• STEP 1. Sample preparation and lysis**

Collect samples in a Collection Plate (first collection plate) → Add VNE Buffer → Mix by pipetting → Stand at room temperature for 5 min

**• STEP 2. Adjust binding condition:**

Add AD Buffer → Mix by pipetting

↓

Centrifuge protocol ✓ or Vacuum protocol

**• STEP 3. Bind DNA/RNA to Filter Plate:**

centrifuge protocol	vacuum protocol
<ul style="list-style-type: none"> <li>• Combine the plates.</li> <li>• Transfer the sample mixture to Filter plate.</li> <li>• Centrifuge at 4,500 – 6,000 x g for 2 min.</li> </ul>	<ul style="list-style-type: none"> <li>• Fix Plates to manifold.</li> <li>• Transfer the sample mixture to Filter plate.</li> <li>• Apply 10 inches Hg vacuum until the wells have emptied.</li> </ul>
<p><b>• STEP 4. Wash the Filter Plate with Wash Buffer 1</b></p> <ul style="list-style-type: none"> <li>• Add Wash Buffer 1. Centrifuge at 4,500 – 6,000 x g for 2 min.</li> </ul>	<ul style="list-style-type: none"> <li>• Add Wash Buffer 1. Apply vacuum at 10 inches Hg.</li> </ul>
<p><b>• STEP 5 &amp; 6. Wash the Filter Plate with Wash Buffer 2</b></p> <ul style="list-style-type: none"> <li>• STEP 5: Add Wash Buffer 2. Centrifuge at 5,600 - 6,000 x g for 2 min</li> <li>• STEP 6: Add Wash Buffer 2. Centrifuge at 5,600 - 6,000 x g for <b>15 min</b></li> </ul>	<ul style="list-style-type: none"> <li>• STEP 5: Add Wash Buffer 2. Apply vacuum at 10 inches Hg.</li> <li>• STEP 6: Add Wash Buffer 2. Apply vacuum at 10 inches Hg for <b>10 min</b>.</li> </ul>
<p><b>• STEP 7. Dry the membranes of the Filter Plate:</b></p> <ul style="list-style-type: none"> <li>• Stand the Filter plate on a clean paper towel at room temperature for 10 min.</li> </ul>	<ul style="list-style-type: none"> <li>• Tap the Filter Plate tips on paper towel</li> <li>• Return the Filter Plate and the Collection Plate to the manifold.</li> <li>• Apply maximum vacuum for an additional 10 min.</li> </ul>
<p><b>• STEP 8. RNA Elution:</b></p> <ul style="list-style-type: none"> <li>• Add RNase-free Water to the Filter Plate. Stand for 3 min.</li> <li>• Centrifuge to elute RNA.</li> </ul>	<ul style="list-style-type: none"> <li>• Add RNase-free Water to the Filter Plate. Stand for 3 min.</li> <li>• Close the manifold valve. Turn on the vacuum source to build up a vacuum to 15 inches Hg.</li> <li>• Open the manifold valve to apply vacuum to elute DNA/ RNA.</li> </ul> <p><b>Alternative: If the consistent volume of elutes are recommend use centrifuge protocol to proceed this elution step. (Page 3, STEP 8)</b></p>

## Protocol: (centrifugation processing)

Please Read Important Notes Before Starting The Following Steps.

### Required hardware

Centrifuge equipment capable of 5,600 ~ 6,000 X g with a swing -bucket rotor and the adaptor for 96-well plate

### STEP 1. Sample preparation and lysis

- Transfer 200 µl of sample to each well of a Collection Plate (provided, 96-well 2 ml plate; first collection plate). If prepared samples are less than 200 µl, adjust sample volume to 200 µl with PBS (not provided).
- Add 400 µl of VNE Buffer to each well and mix completely by pipetting.
- Incubate at room temperature for 10 min.

### STEP 2. Adjust binding condition

- Add 300 µl of AD Buffer (ethanol added) to each well and mix completely by pipetting.

### STEP 3. DNA/ RNA Binding

- Place a Filter Plate (provided, 96-Well nucleic acid binding plate) on a clean Collection Plate (provided, second collection plate).
- Transfer the sample mixture to each well of the Filter Plate and discard the Collection Plate (first collection plate).
- Place the plates in a rotor bucket and centrifuge at 5,600 – 6,000 x g for 2 min.
- Discard the flow-through and return the Filter Plate to the Collection Plate.

### STEP 4. Wash the Filter Plate with Wash Buffer 1

- Add 500 µl of Wash Buffer 1 (ethanol added) to each well of the Filter Plate.
- Place the combined plate in a rotor bucket and centrifuge at 5,600 – 6,000 x g for 2 min.
- Discard the flow-through and return the Filter Plate to the Collection Plate.

### STEP 5. Wash the Filter Plate with Wash Buffer 2

- Add 500 µl of Wash Buffer 2 (ethanol added) to each well of the Filter Plate.
- Place the combined plate in a rotor bucket and centrifuge at 5,600 – 6,000 x g for 2 min.
- Discard the flow-through and return the Filter Plate to the Collection Plate.

### STEP 6. Wash the Filter Plate again with Wash Buffer 2

- Add 500 µl of Wash Buffer 2 (ethanol added) to each well of the Filter Plate.
- Centrifuge at 5,600 – 6,000 x g for **15 min**.
- Discard the flow-through and the Collection Plate (second plate).

### STEP 7. Dry the membranes of Filter Plate

- Place the Filter Plate on top of a clean paper towel (not provided) and stand at room temperature for 10 min.

### STEP 8. DNA/ RNA Elution

- Place a Elution Plate (provided, 96-Well PCR plate) on top of a clean Collection Plate (provided, third collection plate) then place the Filter Plate on the Elution plate. (top: Filter Plate, middle: 96-well PCR Plate, bottom: Collection Plate)
- Add 50 ~ 75 µl of RNase-free Water to the membrane center of the Filter Plate. Stand for 3 min.
- **Important Step! For effective elution, make sure that RNase-free water is dispensed on the membrane center and is absorbed completely.**
- **Important : Do not elute the DNA/ RNA using RNase-free water less than suggested volume (< 50 µl). It will lower the DNA/ RNA yield.**
- Place the plates in a rotor bucket and centrifuge at 5,600 – 6,000 x g for 5 min to elute DNA/ RNA.
- Seal the Adhesive Film and store the RNA at -70 °C.

## Protocol: (vacuum processing)

Please Read Important Notes Before Starting The Following Steps.

### Required hardware

Vacuum manifold for 96-well plate and vacuum source reached to -15 inches Hg

**Alternative:** If using centrifugation for Elution Step (STEP 8), a centrifuge equipment is required, capable of 5,600 ~ 6,000 X g, with a swing -bucket rotor and the adaptor for 96-well plate.

### STEP 1. Sample preparation and lysis

- Transfer 200 µl of sample to each well of a Collection Plate (provided, 96-well 2 ml plate; first collection plate). If prepared samples are less than 200 µl, adjust sample volume to 200 µl with PBS (not provided).
- Add 400 µl of VNE Buffer to each well and mix completely by pipetting.
- Incubate at room temperature for 10 min.

### STEP 2. Adjust binding condition

- Add 300 µl of AD Buffer (ethanol added) to each well and mix completely by pipetting.

### STEP 3. DNA/ RNA Binding

- Fix a clean Collection Plate (provided, second collection plate) on the rack of vacuum manifold and cover the manifold lid. Place a Filter Plate (provided, 96-Well nucleic acid binding plate) on top of the Collection Plate (second collection plate).
- Transfer the sample mixture to the Filter Plate and discard the Collection Plate (first collection plate).
- Apply vacuum at 10 inches Hg until the wells have emptied.
- Discard the flow-through and return the Filter Plate and the Collection Plate to the manifold.

### STEP 4. Wash the Filter Plate with Wash Buffer 1

- Add 500 µl of Wash Buffer 1 (ethanol added) to each well of the Filter Plate.
- Apply vacuum at 10 inches Hg until the wells have emptied.
- Discard the flow-through and return the Filter Plate and the Collection Plate to the manifold.

### STEP 5. Wash the Filter Plate with Wash Buffer 2

- Add 500 µl of Wash Buffer 2 (ethanol added) to each well of the Filter Plate.
- Apply vacuum at 10 inches Hg until the wells have emptied.
- Discard the flow-through and return the Filter Plate and the Collection Plate to the manifold.

### STEP 6. Wash the Filter Plate again with Wash Buffer 2

- Add 500 µl of Wash Buffer 2 (ethanol added) to each well of the Filter Plate.
- Apply vacuum at 10 inches Hg for **10 min**.
- Discard the flow-through and return the Collection Plate to the manifold.

### STEP 7. Dry the membranes of Filter Plate

- Gently tap the tips of the Filter Plate on a clean paper towel to remove residual liquid.
- Return the Filter Plate to the Collection Plate fixed in the manifold.
- Apply vacuum for an addition 10 min.
- Discard the flow-through and the Collection Plate (second plate).

### STEP 8. DNA/ RNA Elution

- Place a Elution Plate (provided, 96-Well PCR plate) on top of a clean Collection Plate (provided, third collection plate) and fix plates on the rack of manifold. Cover the manifold lid and place the Filter Plate on the Elution Plate. (top: Filter Plate, middle: 96-well PCR Plate, bottom: Collection Plate)
- Add 50 ~ 75 µl of RNase-free Water to the membrane center of the Filter Plate. Stand for 3 min.
- **Important Step! For effective elution, make sure that RNase-free water is dispensed on the membrane center and is absorbed completely.**
- **Important : Do not elute the DNA/ RNA using RNase-free water less than suggested volume (< 50 µl). It will lower the DNA/ RNA yield.**
- Close the manifold valve. Turn on the vacuum source to build up a vacuum to 15 inches Hg.
- Open the manifold valve to apply vacuum to elute DNA/ RNA.
- Seal the Adhesive Film and store the RNA at -70 °C.

**Alternative:** If the consistent volume of elutes are recommend use centrifuge protocol to proceed this elution step. (Page 3, STEP 8)