

NEW LYOPHILISED KITS FOR HEPATITIS B AND C WORK SIMPLER AND FASTER!

- Quantitative analysis of samples extracted from blood plasma
- Independent internal control to monitor all test stages (DNA/RNA extraction and amplification) included
- Ready-to-use lyophilized reaction mixture in universal 0,2 ml tubes for quicker and easier preparation
- Room temperature setup
- The same amplification profiles (can be performed simultaneously) of all hepatitis kits (incl. genotyping kits)













HBV

Detection limit 3 IU/ml in a sample

HCV

Detection limit 7 IU/ml in a sample

Catalog number	Description	Detection	Format	# of reactions	Certification	Assay detail	Shelf life	Fluorescent channels	Availability
H-4022-1-14-CE	AmpliSens® <i>HBV</i> -Monitor-L PCR kit variant FRT-L			96		quantitative	12 months		
H-4002-1-14-CE	AmpliSens® <i>HCV</i> -Monitor-L PCR kit variant FRT-L			96		quantitative	12 months		

FINISH THE DIAGNOSIS WITH OUR CE IVD GENOTYPING KITS














HCV GENOTYPES:

1a, 1b, 2, 3a, 4, 5a, and 6

HBV GENOTYPES:

A, B, C and D

Catalog number	Description	Detection	Format	# of reactions	Certification	Assay detail	Shelf life	Fluorescent channels	Availability
H-1982-1-3-CE	AmpliSens® <i>HBV</i> -genotype-FRT PCR kit			55		genotyping	12 months		
H-0256-1-3-CE	AmpliSens® <i>HCV</i> -genotype-FRT PCR kit (REVERTA-L included)	 		55		genotyping	12 months		

-  Real-Time
  Reverse transcription
  Fluorescent End-Point
  Reagents in stock tubes
  Ready-to-use PCR tubes
-  on request
  available

Ver.03_07/2022

